

Cavislim advance

Tecnovita
by **BH**

ENGLISH

FRENCH



cavitation/cavitation

safe and effective
figure shaping

modelage
sûr et efficace
de la silhouette

User
manual
Manuel de
l'utilisateur
YS70

contents

chapter 1

- contents of the box
- assembling CAVISLIM ADVANCE

chapter 2

- basic characteristics of CAVISLIM ADVANCE
- operating instructions

chapter 3

- application modes

chapter 4

- correct, safe use of CAVISLIM ADVANCE

chapter 5

- frequently asked questions about CAVISLIM ADVANCE

chapter 6

- important safety instructions
- warranty
- certificates
- technical specifications
- parts list


chapter 1

EN


contents of the box

Inside this box you will find everything you need to enjoy CAVISLIM ADVANCE from now onward. Carefully check that the pack contains all of the elements on the following list .

Hand unit support



Hand unit M80




Hand unit M38



Main unit




Operating pedal



DVD



User manual



Cavitation gel




Allen key



4 screws



Mains cable



Anti-cellulite serum




Toning serum



assembling Cavislim Advance

First we are going to fit the **handset support** onto the **main unit**. To do this, place the main unit on a table and position the **handset support** by inserting the **screws** as shown in the diagram.

Using the **Allen key** tighten the **screws** securely so that the two parts fit together perfectly.



chapter 2

EN

basic characteristics of Cavislim Advance

CAVISLIM ADVANCE is a hi-tech device that allows us to selectively and efficiently remove the build-up of fat (cellulite, orange skin syndrome) from our body. The technology used is that of ultrasound emission which produces the physical phenomenon known as cavitation, causing the removal of fat cells.

This technology with proven results is widely used in beauty salons all over the world in order to fight cellulite painlessly and non-invasively. Characteristics for model CAVISLIM ADVANCE Ref. YS70

The YS70 model includes 2 handsets for performing 2 different types of treatment.

The larger **M80** handset is for general, all-over body treatment. It has 4 pulsed mode programs and 1 continuous mode program, with 5 different intensity levels.

Time: Programmable from 1 to 15 minutes.

Intensity: 5 intensity levels (5-10-15-20-25)

Mode: 4 programmable pulsed work modes (P1-P2-P3-P4)
1 continuous work mode (P5)

The smaller **M38** handset is specially designed to treat the more troublesome areas of fat, because the treatment area is smaller the power is more concentrated and its effect is greater. It has 10 pulsed mode programs and 4 different intensity levels. There is also an option to light an LED ring for greater control over the area being treated.

Time: Programmable from 1 to 15 minutes.

Intensity: 4 intensity levels (1-2-3-4)

Mode: 10 programmable pulsed work modes (P1-P2-P3-P4-P5-P6-P7-P8-P9-P10)

With this unit you only have to worry about programming in three simple functions: intensity, time and work mode.

The intensity function controls how deep the sound waves penetrate your skin. The greater the intensity, the deeper the waves reach.

The time function allows you to choose how long you wish to work with CAVISLIM ADVANCE. For reasons of safety the time function allows you to choose up to 15 minutes which corresponds to the maximum recommended time for working a specific area. However, you can reduce this maximum time by 1 minute intervals. Once the selected time has elapsed CAVISLIM ADVANCE will alert you with a beep and a message on the display showing the text "Time out".


CAVISLIM ADVANCE also lets you choose the work mode, helping you select the time interval during which the waves reach your skin.

Cavislim Advance operating instructions

Precautions: Consult your doctor before beginning to use this appliance. This warning is especially important for people with health problems. Use by those under 18 years of age is not recommended. Keep this appliance out of children's reach.

Starting


- Connect the **handsets** to the main unit. To do this, take the plug for the **M38 handset** and insert into the socket on the left side of the **main unit** as shown in Fig. 2.1
- Next, take the plug for the **M80 handset** and insert into the socket on the right side of the **main unit** as shown in Fig.2.2
- Now take the jack for the **operating pedal** and insert into the socket on the left side of the main unit as shown in Fig. 2.3 Check that all the connections are made correctly before connecting it to the mains supply.
- Finally, plug one end of the **mains cable** into the socket on the machine and the other end into the electric wall socket (220V) Fig. 2.4 Switch the On/Off switch to position (I) to start the machine.



Monitor elements

The CAVISLIM ADVANCE monitor consists of a display screen, an LED and five buttons or keys.


- The display provides information about the cavitation session.
- The ON LED (green) tells you that the device is operating.
- The keys have the following functions:
 - ▲ and ▼ change the work **MODE**
 - + and - adjust **TIME** and **INTENSITY**
 - || pause and ► start session



Cavitation session

There are 3 settings for the massage:

- **Time:** measures the duration of the session in minutes and seconds. You can select between 1 to 15 minutes (in steps of 1 minute).
- **Work mode:** Pulse type (impulse emission frequency). CAVISLIM ADVANCE allows you to choose different work modes depending on the handset being used:
 - **Large handset** (CH = M80) offers a choice of 5 modes: P1, P2, P3, P4 and P5 (P1 very low – P4 very high, P5 continuous)
 - **Small handset** (CH. = M38) offers a choice of 10 modes: P1, P2, P3, P4, P5, P6, P7, P8, P9 and P10.
- **Intensity:** establishes the intensity level for the session. CAVISLIM ADVANCE allows you to choose different intensity levels depending on the handset being used:
 - **Large handset** (CH. = M80) offers a choice of 5 intensity levels: 5, 10, 15, 20 and 25.
 - **Small handset** (CH. = M38) offers a choice of 4 intensity levels: 1, 2, 3, 4.



Beginning a session

Once you have plugged the machine into the mains supply and switched it to position (I) the following message will appear on the monitor: **WELCOME**

This message will remain visible for 3 seconds, after which the following default settings will appear:

TIME = 10:00
CH. = M80
MODE = P1
INT. = 5

From here you can modify the time, handset and work mode settings to establish the type of session you want:

Select the parameter to be modified with the small arrow located initially next to the TIME concept which you can move it with **+** and **-** keys

Changing session duration

Once the TIME concept is selected, you can modify the session duration with the **+** and **-** keys.


The displayed default value is 10 minutes but this can be changed in steps of 1 minute between the values 1:00 and 15:00.


Modify the time parameter at any moment during the session.

Selection of handset.

Once the CH concept is selected, the handset to be used can be changed in the session with the **+** and **-** keys.

If you choose the M38 handset you can select if you want to switch on the red lights on its head or not, by pointing with the small arrow concept **P**. With the **+** and **-** keys, you can choose if the lights should be switch on, by letter Y, or if it should be permanently off, by letter N.






Changing work mode

Once the MODE concept is selected, you can modify the work mode with the **+** and **-** keys.

The default value shown is P1 (very low pulse frequency) but this can be modified by selecting from the values **P1**, **P2**, **P3**, **P4** and **P5** (with handset **M80**) and **P1**, **P2**, **P3**, **P4**, **P5**, **P6**, **P7**, **P8**, **P9** and **P10** (with handset **M38**).



Changing message intensity

Once the INT. concept is selected, you can modify the message intensity with the **+** and **-** keys.

The default intensity level is 5 but this can be changed by selecting from levels 5, 10, 15, 20 and 25 (with handset **M80**) and levels 1, 20, 3 and 4 (with handset **M38**).

Finally, press the ► button and with the machine on (a green light appears in the centre of the unit), take hold of the handset, place it on the part where you want to work and press down on the pedal to begin your cavitation session.

When you tread on the pedal you will see the time count down on the screen and the handset will begin to emit impulses. Place the handset on the area you wish to work and massage it as shown on the DVD on in the “Application modes” section.


If you want to pause the session, just stop pressing on the pedal and the impulse emission will stop immediately.

When the time counter reaches **0.00** the unit will stop and emit a series of beeps to tell you that the selected session has finished. The following message may appear on the monitor:

**TIME UP
PRESS A KEY
TO START**

Pressing any key on the machine will take you back to the default screen.

**TIME = 10:00
CH. = M80
MODE = P1
INT. = 5**



chapter 3

application modes

Always apply an abundant amount of gel to the area that you wish to work and begin to treat it by sliding the handset slowly and without pressing, as shown in the diagrams.


If you need to purchase gel or any of our preparation creams you can find them in the shop where you purchased CAVISLIM ADVANCE or at our website: www.tecnovita.eu

In each of the following sections you will find a table with application protocols. Before beginning to use the appliance, you must determine the level of fat in each area to be treated.

- A** Large accumulation of high density fat
- B** Large accumulation of medium density fat
- C** Medium accumulation of medium density fat
- D** Small accumulation of fat
- E** orange peel syndrome
- F** Maintenance session

central abdomen

Apply an abundant amount of gel on your abdomen and begin to massage the area by sliding the handset slowly and without pressing, as shown in the diagrams.




For best results we recommend applying CAVISLIM ADVANCE by combining different massage techniques (spiral, crossways, circular) in different sessions. We recommend lying down for a more comfortable CAVISLIM ADVANCE application.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	3	15	9	3	10'
B	3	15	9	3	10'
C	3	10	8	2	7'
D	2	10	8	1	7'
E	2	10	7	1	7'
F	2	10	6	1	7'

Remember, never apply CAVISLIM ADVANCE to the navel, breasts, the thorax area or the genitals.

sides

EN




Apply gel over the area to be worked using CAVISLIM ADVANCE

Circular and longitudinal massage movements should be applied from top to bottom as shown in the diagrams. It is important to massage in the correct direction, as this will improve drainage. Combine the various massage techniques during different sessions.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	4	20	7	3	10'
B	3	15	7	3	10'
C	2	10	6	2	7'
D	1	5	4	1	7'
E	1	5	4	1	5'
F	1	5	4	1	5'

hips



Apply gel over the area to be worked using CAVISLIM ADVANCE.


Circular and longitudinal massage movements should be applied from the bottom upward as shown in the diagrams.

During treatment, concentrate on the area with the most fat.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	3	15	9	3	7'
B	2	15	8	3	7'
C	2	10	8	2	7'
D	1	5	6	2	5'
E	1	5	4	1	5'
F	1	5	4	1	5'

buttocks

Apply gel over the area to be worked using CAVISLIM ADVANCE.




Both longitudinal and circular massage should be done in the direction of the arrows shown in the diagrams.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	5	25	9	4	10'
B	5	20	8	3	7'
C	4	20	7	3	7'
D	3	15	6	2	7'
E	2	15	6	2	10'
F	2	15	5	2	10'

inner thigh

Remember, never apply CAVISLIM ADVANCE from the knee to the foot, there is no significant accumulation of fat in this area. Apply gel over the area to be worked using CAVISLIM ADVANCE.



Both the longitudinal and circular massage should be done from the bottom upward as shown in the diagrams. Follow these simple instructions to help you get the best results.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	3	15	8	2	10'
B	3	15	8	2	10'
C	2	10	6	2	7'
D	2	10	6	2	7'
E	2	10	4	2	7'
F	2	5	3	1	7'

inner knee

EN




Apply a light circular massage to this area, as shown in the diagram.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	3	5	7	2	10'
B	3	5	7	2	10'
C	2	5	6	1	7'
D	2	5	4	1	7'
E	1	5	1	1	7'
F	1	5	1	1	7'

Recuerde no aplicar el CAVISLIM ADVANCE sobre las articulaciones de la rodilla.
Aplique gel sobre la zona a trabajar con el CAVISLIM ADVANCE.

Remember not to apply CAVISILIM on the knee joints.
Apply gel over the area to be worked using CAVISLIM

arms



Put a generous amount of gel on the area to be worked and this will improve the effectiveness of the treatment and prevent damaging the machine. Perform the massage, sliding the handset slowly and without pressing. Massaging this area should be applied toward the body, as shown in the diagram.

Nivel	M80		M38		Time
	Mode	Intensity	Mode	Intensity	
A	2	5	6	1	7'
B	2	5	4	1	5'
C	1	5	2	1	5'
D	1	5	2	1	3'
E	1	5	1	1	3'
F	1	5	1	1	5'

Recuerde nunca aplicar CAVISLIM ADVANCE en la zona de manos , antebrazo y articulaciones

Remember, never apply CAVISLIM ADVANCE on the

chapter 4

correct, safe use of Cavislim Advance

For correct and safe use of CAVISLIM ADVANCE, you should keep the following application times in mind:

- **Per work area: 15 minutes maximum.**

By "work area" we mean each one of the body parts where CAVISLIM ADVANCE can be applied. The various work areas are: arms, abdomen, waist, hips, inner thighs and inner knees.

E.g.: In a single session you can apply CAVISLIM ADVANCE for 15 minutes on your sides and 15 minutes on your arms. But, for example, you must NEVER perform a 30 minute application session just on your sides in the same day.

- **By session we mean using CAVISLIM ADVANCE on two different work areas (15 minutes per area) in the same day.**

- **Between sessions: 72 hours minimum.**

By "between sessions" we mean the rest time between each session. E.g.: In one week the session schedule (maximum) would be: Monday, Thursday and Sunday.

- **Every 20 sessions:**

You should rest for 3 weeks. After this period you may decide to repeat the first treatment or start the maintenance program.

- **Maintenance sessions:**

This depends on the individual, but one session per month should suffice.

VERY IMPORTANT

- **Never exceed the application times (15 minutes per area, 30 minutes total in one day).** Your body does not know "where" you are eliminating the fat from, so a session must never include 15 minutes on your sides, 15 minutes on your legs, 15 minutes on your arms, etc. REMEMBER: the maximum time during one day, adding all the parts of the body being treated, must be no more than 30 minutes.

- **And, above all, respect the rest period of 72 hours between sessions.** This is the time that your body needs to eliminate all the fat dissolved during the session through the lymphatic system.

chapter 5

EN

frequently asked questions

What is cavitation?

Cavitation is a non-invasive (surgery free) body shaping technique that reduces volume by eliminating cellulite and localised fat deposits in the abdomen, hips, flanks, arms, inner thighs and knees.

What are the advantages over other methods?

The main advantage is that we can choose the specific areas of the body that we want to work. It could be said that it is a technology which allows us to concentrate on cellulite selectively.

Is cavitation painful?

Use of CAVISLIM ADVANCE should never cause pain.

Is it normal to hear a slight buzzing in your ears while using CAVISLIM ADVANCE?

Yes, some users may detect a slight buzzing sound during its application.

Can it be used on all parts of the body?

It cannot be applied to all parts of the body. Under no circumstances can it be applied to: the head or face, neck, breasts, the thorax, navel, forearm, joints, genitals, back, ankles or lymph glands.

Does CAVISLIM ADVANCE reduce weight and volume?

CAVISLIM ADVANCE eliminates clusters of fat, which have weight; consequently your body loses volume.

Is CAVISLIM ADVANCE effective on orange peel syndrome?

Yes, CAVISLIM ADVANCE is effective on orange peel syndrome because it destroys adipose cells that cause fat storage.

Can CAVISLIM ADVANCE be used to remove wrinkles?

No, CAVISLIM ADVANCE has not been designed to remove wrinkles.

Can CAVISLIM ADVANCE be used on the face?

CAVISLIM ADVANCE must NEVER be used on the face, eyes or neck.

After how many sessions will I begin to see results?

You should notice results after the first five sessions, however results may vary depending on the subject.

How long should I wait between sessions?

It is advisable to wait 72 hours between sessions.

Does cellulite reappear after using CAVISLIM ADVANCE?

Although CAVISLIM ADVANCE ensures effective and long-lasting results, we cannot guarantee against cellulite reappearing. To prevent the appearance of new cellulite we recommend a balanced diet, regular

Are there any side effects?

If you are healthy, CAVISLIM ADVANCE has no reported side effects; nevertheless in some people it may produce a reddening of the skin that can last for about 2 hours.

Is it necessary to combine CAVISLIM ADVANCE with a special diet?

No, it is not necessary to follow any special diet, but a balanced diet and proper hydration is highly advisable.

How much water should I drink?

It is essential to begin hydrating yourself two hours before applying CAVISLIM ADVANCE. It is also advisable to drink more water after finishing a session. The recommended quantity is at least 2 litres during these periods, as it helps to eliminate toxins.

Is it necessary to combine CAVISLIM ADVANCE with a sports activity?

No, this is not strictly necessary, but at least 30 minutes of physical activity will enhance the results as it helps to drain fluids.

Which sports activities can supplement CAVISLIM ADVANCE?

The vibration platform is the ideal complement due to its toning and draining qualities. But if you do not have access to a platform, then walking, exercise bike, swimming or any other aerobic exercise for at least 30 minutes is very helpful.

How can I enhance the effects of CAVISLIM ADVANCE?

By combining a balanced diet with regular physical exercise and drinking plenty of water. Added to which, drainage massage can also be very useful.

What are the anti-cellulite and toning gels for? When do I use them?

1. Anti-cellulite gel: Contains active ingredients that help eliminate cellulite.
2. Toning gel: Contains active ingredients that make the treated area firmer.

The application sequence is as follows: during the first few days of treatment use the anti-cellulite gel and once you begin to notice results, use the toning gel. This is applied before the Cavitation gel, as with any other serum. The area should be massaged until it is well absorbed and then apply plenty of cavitation gel to start the session.

Can I use CAVISLIM ADVANCE without applying gel?

NEVER, this is because the gel acts as a conductor for the ultrasonic sound waves with your body. If you do not use gel then the CAVISLIM ADVANCE application loses its effectiveness and the machine can be severely damaged.

Can I use any type of gel?

We recommend using only TECNOVITA by BH gel, as the use of any other type may damage the machine.

When is it inadvisable to use CAVISLIM ADVANCE?

- If you are pregnant or think you might be.
- If you are breast-feeding.
- Undergoing fertility treatments.
- If you have any kind of heart disease.
- If you suffer from skin problems.
- Osteoporosis.
- Fibromyalgia.
- Malignant tumours.
- Hypertension.
- High cholesterol levels.
- Kidney problems
- Physical shortcomings
- Metal implants or prostheses.
- Physical malformations.
- Epilepsy.
- Long-term treatment with medications: anti-inflammatories, anticoagulants, antibiotics or antihistamines.
- Depression.
- Psychiatric or psychological treatment.
- Liver problems (fatty liver)
- Irritable bowel syndrome.
- Diabetes

Can I use CAVISLIM ADVANCE if I am pregnant?

You must never use CAVISLIM ADVANCE if you are pregnant or think you might be.

Can CAVISLIM ADVANCE be used during lactation?

Do not use it during lactation.

Can CAVISLIM ADVANCE be used on the breasts?

NEVER apply it to the breasts.

Can I use CAVISLIM ADVANCE during menstruation?

Use on the abdominal area is not advisable during menstruation as it may cause increased bleeding. However, it can be applied in other areas.

Can I use CAVISLIM ADVANCE if I have an IUD?

Yes.

Can CAVISLIM ADVANCE affect my cholesterol levels?

Use of CAVISLIM ADVANCE may increase your cholesterol levels.

Can CAVISLIM ADVANCE be used with skin conditions?

It should not be used if you have any skin condition. Consult your doctor.

Can I use CAVISLIM ADVANCE if I have ear implants?

Do not use it if you have implants fitted.

Can I use CAVISLIM ADVANCE when I am bathing or showering?

Under no circumstances should the handset or the machine come into contact with water due to the risk of electrocution.

How do I clean CAVISLIM ADVANCE and the handset?

Never immerse CAVISLIM ADVANCE or the handset in water. Use a cloth slightly dampened with water.

chapter 6

important safety advice

This product is supplied with a mains cable that has an earth wire and a plug with an earth pin.

In the event that the appliance malfunctions then the earth connection will redirect the electrical current away from the machine and reduce the risk of an electrical shock.

WARNING! If the equipment is not earthed correctly then there is a risk of electrical shock.

Do not modify the plug supplied with the equipment, if it does not fit the wall socket then get an electrician to install an electrical socket that does. This product should be used on a 220-240 volt circuit and has an earthed plug similar to that shown in Fig. 6.1

Follow these recommendations:

1. Never use accessories that are not recommended by the manufacturer.
2. Do not use the machine if the mains cable becomes damaged or worn.
3. Keep the mains cable away from hot surfaces.
4. Do not use the machine outdoors.
5. To unplug the machine, first put the switch to the «0» position and then unplug it from the wall socket.



PROTECTION

This model is equipped with an electronic circuit breaker to protect both the electronic and electrical circuits.

This safety device (fuse) will be automatically triggered if the unit is subjected to abnormal force or heat.

Should this occur:

1. Unplug the cable from the mains supply.
2. Set switch (L), located on the casing, to the «0» (OFF) position. Fig. 6.2
3. Change the fuse (C) (1 Amp).
4. Plug the unit into a wall socket equipped with an EARTH terminal. Set switch (L) to position [1]. Fig. 6.2



warranty

EN

We suggest that you read through this document carefully in order to fully understand the cover that the manufacturer offers as well as your responsibilities as owner.

GUARANTEES

GENERAL INFORMATION

EXERCYCLE, S.L. would like to thank you for purchasing this product and the trust you have placed in us. Our products have been subjected to rigorous quality checks. Nevertheless, should you experience any problem with this appliance during the guarantee period then please read through these warranty conditions and then get in touch with your local authorised dealer. This guarantee certificate has been drafted in accordance with that stated in Directive 1999/44/EC of the European Parliament and of the Council 25th May, 1999, regarding on certain aspects of the sale of consumer goods and associated guarantees, incorporated into Spanish Law by virtue of Law 23/2003, 10th July; Guarantees on the sale of consumer goods.

To this effect, EXERCYCLE, S.L. VAT nºbc B-01029032, with business address c/ Zurrupitieta 22, Vitoria, Spain, guarantees this product under the following terms and conditions:

1. GUARANTEE PERIOD

The product covered by this guarantee certificate is protected for a period of TWO (2) YEARS as from the date of receiving the product, which will be taken as the date that appears on the sales receipt for the product or on the corresponding delivery docket should this be later. It is essential that this be presented in the event of possible repairs. The policy regarding the warranty period may vary between different countries. The stated warranty period is only applicable to Europe.

2. RIGHTS UNDER THE GUARANTEE

- Repair or replacement of the product which, in such instance, will be carried out within a reasonable period and without significant inconvenience to the consumer, taking into account of the nature of the goods and the purpose for which the consumer required them.
- Those cases in which the repair or replacement represents a disproportionate cost for the seller in comparison with other forms of remedy will be exempt assuming these do not prejudice the consumer.
- Repairs will be free for the consumer with regard to labour, materials and carriage.
- In the event that the machine is repaired or replaced with a new one, the guarantee period will be suspended temporarily until the new appliance has been delivered and will carry at least six months of guarantee, however under no circumstances will a new 2 year period be initiated.
- The buyer may, as is apparent, require (i) a reduction in price or proceed to a (ii) recession of the contract when the repair or replacement is not carried out within a reasonable time or with significant inconvenience to the buyer.

The reduction in price (i) will be proportional to the difference between the value that the goods would have had if they had been delivered in accordance with the contract and the actual value that the goods were worth at the time of the said delivery. Recession of the contract (ii) will not proceed when the non-conformity is of minor importance.

3. THE GUARANTEE DOES NOT COVER

- Defects attributable to use other than that established in the User Manual.
- Professional use or commercial leasing or hiring of the product.
- Poor installation, poor assembly or breakage of components due to improper installation, (the User Manual gives detailed instructions on how to install and use the machine).
- Improper use of the product (use without applying gel, exceeding operating times, exposure to high moisture levels, etc.) or problems arising from poor maintenance (accumulation of dirt on any part of the product or its accessories, etc.) (check the User Manual for detailed instructions on how to install and use the product).
- Repairs carried out by persons not authorised by EXERCYCLE, S.L. or appliances that have been opened or tampered with by persons not authorised by the manufacturer.
- Deviations from the standard product, such as light vibration, noise or breakage of components that do not affect the performance of the machine and which cannot be attributed to manufacturing faults.
- Scratches, scrapes, abrasions, knocks, falls, etc. caused by carelessness.
- Defects on the appliance caused by damage during transit.

4. CLAIMING AGAINST THE GUARANTEE

- Claims made for nonconformities will be limited to three years as from the date that the goods were delivered.
- The consumer must report the existence of the defect within a period of 2 months as from the date that the defect was discovered.
- This guarantee does not affect the consumer's statutory rights under any applicable law, unless it is incompatible with the application of other actions derived from compensation sought for hidden faults in the contract of sale.
- With regard to queries that may arise from the interpretation or application of this guarantee, these will be governed by the corresponding national law.


Any discrepancy arising from the interpretation, validity or application of this contract will be subject to the courts and tribunals of Vitoria, Spain.

certificates

By virtue of this document, Exercycle S.L. declares that this product meets the essential requirements and other relevant provisions of Directives 2004/108/EC and 2006/95/EC

parts list

1 Operating
pedal




2 Handset
M80




17 Handset
M38




3 Handset
support



G01 Main unit



16 Mains
cable



To order replacement parts, please indicate:
The machine model / Corresponding part N°/ Quantity

N°	Código
1	YS50001
2	YS50002
3	YS70003
16	YS50016
17	YS70017
G01	YS70G01

Example:

YS50001	1
---------	---

technical specifications

Net weight: 7 kg

Dimensions: 34 cm L x 38 cm W x 13 cm H

Ultrasound: 40 kHz

Output power: 25 Watts maximum (depending on program)

Power supply: Voltage: 220 V AC / 50 Hz frequency